



AED MONTHLY MAINTENANCE CHECKLIST

Inspected by:

Date

- AED is in a clear and visible location with no obstructions. Signs clearly mark location.
- Visual inspection to make sure there is no obvious damage
- Check AED battery by making sure the light is blinking or there is a green check mark.
Battery Expiration Date: _____
- Are the AED pads in good condition and unexpired? Pads Expiration Date: _____
- Do you have both adult and pediatric pads?
- Are gloves, CPR barrier device, scissors, razors, and gauze/ cloth present? (These items should be in your fast response kit.)

Notes:



**First Response
Safety Training**

(770) 945-6837
Info@firstresponsecpr.com
www.firstresponsecpr.com



THINGS TO-DO

Name

Date

To-do

- _____
- _____
- _____
- _____
- _____
- _____
- _____
- _____
- _____
- _____
- _____
- _____

Notes

Goal